Colourflow

lowable dental composite

MD Medical device class IIa

COLOURFLOW is a flowable polymer dental material type 1 (class 2, (aluminium has a radio-opacity equivalent to that of dentine; thus 1 mm group 1), which meets the requirements of the ISO 4049 standard. COLOUR-ELOW composite undergoes free radical polymerisation activated by visible light from the blue region (400-500 nm)

COLOURFLOW has a radio-opacity equivalent to 4-5 mm of aluminium

INDICATIONS FOR USE

- Black's class I, II, III, IV and V cavities
- sealing of pits, fissures and cavities in deciduous and permanent teeth

INSTRUCTIONS FOR USE

FILLING CAVITIES IN DECIDITOLIS AND PERMANENT TEETH

Etch using ETCHGEL, rinse and gently dry the surface (A). Thoroughly cover the surface with a thin layer of a bonding system, blow off the excess and polymerise (B). In the case of very deep cavities cover the cavity bottom with a thin laver of liner and then apply lavers of COLOURFLOW composite (approx, 1 mm thick).

Cure each laver (C) in accordance with the polymerisation table provided. Cure in sections not longer than 5 mm at a time! Adjust to the bite and polish.

- first contrasting composite layer in the case of temporary restorations marking root canal orifices after endodontic treatment
- of material having a radio-opacity equivalent to 1 mm of aluminium has a radio-opacity equivalent to that of dentine and 2 mm of aluminium is equivalent to enamel)





SFALING

Etch using ETCHGEL, rinse and gently dry the surface (A). Thoroughly cover the surface with a thin layer of a bonding system (B), blow off the excess and polymerise

Apply thin layers of COLOURFLOW flowable composite (approx. 1 mm). Cure each layer (C) in accordance with the polymerisation table. Cure in sections not longer than 5 mm at a time!



R

ROOT CANAL MARKING

It is advised to use COLOURFLOW composite to mark different conditions of root canals in order to facilitate diagnosis during subsequent treatment. For the suggested root canal marking system visit our website www.dentallifesciences.com



FIRST CONTRASTING COMPOSITE LAYER IN THE CASE OF TEMPORARY RESTO-RATIONS

TEMPORARY SPLINTS

۸

Etch and apply a bonding system (such as DLS BOND) in accordance with the instructions. Apply a small amount of COLOURFLOW on the teeth intended for splinting and cure each tooth separately in accordance with the polymerisation table. The layer of COLOURFLOW should be clearly visible to enable an easy distinction between the dentine and the composite which is later to be removed together with the fibre used for splinting. Do not use COLOURFLOW black and red as the first composite laver

in the case of temporary splinting.





TEMPORARY OCCLUSAL CORRECTIONS AND MAINTENANCE OF INTERDENTAL SPACES





DLS DENTAL LIFE

POLYMERISATION TABLE

Lamp	COLOUR	Polymerisation depth dependent on duration of the exposure	
		20 s	30 s
Halogen/LED (500-800 mW/cm ²)	pink, green, blue, yellow, purple	1.5 mm	2.0 mm
	white, orange	-	1.0 mm
LED (> 800 mW/cm ²)	pink, green, blue, yellow, purple	2.0 mm	2.5 mm
	white, orange	1.0 mm	1.5 mm
	red, black	-	1.0 mm

Attention! Do not use halogen lamos for the polymerisation of COLOURFLOW red and black! Polymerise COLOURFLOW red and black in thin layers, not thicker than 1 mm.

COMPOSITION

ADVICE FOR DENTAL PERSONNEL

Mixture of dimethacrylate resins: BisGMA, TEGDMA, UDMA, BisEMA: mineral fillers (approx, 61% by weight); Al-Ba-B-Si glass, Ba-Al-B-F-Si glass, pyrogenic silica, pigments; photoinitiation system (CQ : DMAEMA).

Inorganic filler particle size is in the range of 20 nm to 2.0 µm. COLOURFLOW can be used with any standard, light cured, dimethacrylate resin-based bonding systems. A bonding system should be applied in accordance with the manufacturer's instructions.

CONTRAINDICATIONS

Do not use COLOURFLOW composite in patients with a known acrylates allergy

Do not use the product in patients with a hypersensitivities to any of the components.

ADVERSE REACTIONS

None known. However, an allergic reaction cannot be excluded in particularly sensitive individuals.

LIMITATIONS IN USAGE. INTERACTIONS

Do not use with materials containing phenolic compounds, especially eugenol and thymol. Such materials may disrupt polymerisation of the composite. Do not use if it is impossible to completely isolate the area from saliva, blood or moisture. Contamination may disrupt the polymerisation process. Do not use if the syringe or the applicator are suspected to be defective or damaged.

Do not use when any change in product properties is found.

PRECAUTIONS FOR PATIENTS

This device contains substances that may cause an allergic reaction in certain individuals. Do not use in patients with a known acrylates allergy. Avoid contact of an uncured product with skin, eyes and soft tissues of the mouth. If a prolonged contact occurs, rinse with plenty of water. If an allergic reaction occurs, seek medical attention as needed; remove the product if necessary and discontinue future use of the product.

In case of swallowing or aspiration into the respiratory tract, seek immediate medical attention.

If any changes in the work are noticed, attend a dental check-up.

PRECAUTIONS FOR DENTAL PERSONNEL

This device contains substances that may cause an allergic reaction in certain individuals. To avoid the risk of such reaction, minimise the contact with an unpolymerised composite. If contact with skin occurs, rinse with plenty of water. To minimise the risk of contact, always wear personal protective equipment such as gloves, face masks and safety glasses. Acrylates may penetrate some commonly used gloves.

If any contact with a glove occurs, remove the glove and discard it; wash your hands with soap and water and put on a new glove. If an allergic reaction occurs, seek medical attention as needed. The applicators provided with the syringe are blunt in order to reduce the risk of injury, but they should always be handled with care.

To isolate the operative field and to protect the patient the use of a rubber dam is recommended

Ensure sufficient polymerisation of the entire composite layer. Insufficiently, polymerised product can be allergenic and the lifetime of such work may be shortened. In case of insufficient polymerisation, remove the incorrectly cured layer and apply another one, curing it correctly.

In case of any contamination of an uncured composite, the contaminated layer must be removed. In case of contamination or mechanical damage to an already polymerised layer, gently etch its surface and reapply the composite. Cure in accordance with the polymerisation table provided. To minimise the risk of the potential release of unwanted substances, always clean and rinse the surface immediately after curing. Examine the work at a preventive visit or a check-up. In case of any changes in the performance of the work (e.g. wear, chipping), remove the defective work and replace it with a new one. Inform the patient about the need to maintain proper oral hygiene.

WARNINGS

Avoid contamination of the syringe surface (the risk of cross infection). The syringe cannot be reprocessed using heat sterilisation or immersion in a high-level disinfectants. Do not reuse the syringe if it becomes contaminated. If you apply COLOURFLOW directly from a syringe. use one applicator for one patient only due to hygiene reasons.



Ensure the syringe is recapped properly after each use. Polymerisation of the composite may be initiated by ambient light or by a dental operating lamp. To avoid accidental polymerisation of the composite in the applicator, always pull back the syringe plunger immediately after use.

Keep out of reach of children and unauthorised persons. Use in accordance with the manufacturer's instructions.

Do not use after the expiry date.

STORAGE

Protect against mechanical damage. Store at a temperature under 30°C. If stored at a lower temperature, bring back to room temperature before use. Do not expose to direct sunlight. Protect from light. Do not overheat. Do not freeze.

For use by dentists, dental hygienists and dental technicians only.

PACKAGE CONTENTS

1 syringe with a luer-lock cap (1 g) and 1 disposable applicator.

WARRANTY

ARKONA will replace products that have been proved to be defective or will refund the price of purchase. ARKONA is not liable for any loss or damage caused by misuse or improper use of the product.

Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the state in which the user and / or patient is established.



ARKONA

ARKONA Laboratorium Farmakologii Stomatologicznej Nasutow 99C, 21-025 Niemce, Poland, EU www.arkonadent.com

Instruction for use issued on: 18.05.2022

