



easy to open



ready to use



easy to load



safe disposal

Totally eliminates the risk of cross infection



Notes

.....

.....

.....

.....

.....

.....

.....

.....

.....

.....

.....

.....

.....

.....

.....

.....

.....

.....

.....

.....

.....

.....

.....

.....

NEW

UNIQUE



CROSS INFECTION CONTROL
COMPOSITE



ONE PATIENT DOSE

**CONVENIENCE
AT LOWER COST**


DLS DENTAL LIFE SCIENCES

DENTAL LIFE SCIENCES (mfg) Ltd.
Unit A3, Wigan Enterprise Court,
Wigan Enterprise Park, WN2 2LE
Manufactured in the United Kingdom

www.dentallifesciences.com

CROSS INFECTION CONTROL COMPOSITE


Filler content	78 wt%
Flexural strength [MPa]	135 ± 10 MPa
Water sorption	16 ± 2 µg/mm ³
Water solubility	0.7 ± 0.1 µg/mm ³



Exemplary results of the flexural strength analysis performed in accordance with ISO 4049

Depth of cure

Lamp	Shade	20 s	30 s
Halogen/LED (500-800 mW/cm ²)	A1, A2, A3, B1, B2, C2, D2	2.0	2.5
	A3.5, OA2, OA3, D3	1.8	2.2
	Incisal	3.0	3.5
LED (800-1200 mW/cm ²)	A1, A2, A3, B1, B2, C2, D2	2.5	3.0
	A3.5, OA2, OA3, D3	2.2	2.5
	Incisal	>3.5	>3.5



Even distribution of filler particles (SEM image x 5000)

**QUALITY AND PERFORMANCE -
NO COMPROMISE**

Data on file DLS research laboratory.

Tailored to the way you work - a system created with you in mind

Do you work alone or with an assistant?

With an assistant



15 blisters

20 blisters

If you work with an assistant our packs are designed for easy access and clear indication of how many blisters are left to use.

Alone



Vertical windows
let you see
when to replace

10 blisters

Dispenser - your third hand

A2

The shade is clearly marked on each blister.

Shades: A1, A2, A3, A3.5, B1, B2, C2, D2, D3, OA2, OA3, Incisal