

COMCORD - BRIDGE





A single tooth missing (25).



Fillings from the adjacent teeth have been removed and spaces for bridge abutments have been prepared.



COMCORD connecting the outermost points of the preparation. The fibre is placed on the composite used for the reconstruction of the gingival surface.



Gingival papillae secured by interdental wedges which also facilitate the modelling of interdental surfaces of the tooth that is being reconstructed.



The finished bridge (an occlusal surface view).



The finished bridge (a buccal surface view).

D L S DENTAL LIFE
SCIENCES



www.dentallifesciences.com